

MK-H2300

Surgical Headlight



*Better care
Better comfort!*

MK-H2300

SURGICAL HEADLIGHT

DESCRIPTION

The MK-H2300 LED Surgical Headlight, tailored for ENT, Dental, Vet, Plastic surgery, and more, delivers exceptional brightness and enduring battery life. Its headband version, compatible with surgical loupes, ensures comfortable wear for doctors, offering precision in diverse medical applications.

FEATURES

- The latest Cold LED technology makes the illumination better brightness max. 65,000Lux
- Comparable with daylight in high color rendering index ≥ 93
- Clip-on which provides a longer working time of 5-10hrs
- Wireless integrated
- Light Weight for headband 193g
- Significantly less energy consumption and heat generation
- Headband version, can be connected with surgical loupes and provides comfortable wear for doctors

| Model | MK-H2300 |
|--|---|
| Luminance (LUX) | Up To 65,000 |
| Illumination Spot diameter (D= Φ 200MM) | 10-45mm |
| Illumination Spot diameter (D= Φ 300MM) | 10-75mm |
| Illumination Spot diameter (D= Φ 500MM) | 20-130mm |
| Color temperature (K) | 4300-4800 |
| LED life | 50000hrs |
| Battery Life | 5-12 hrs |
| Color Rendering Index | >93 |
| Magnification Glasses Optional | 2.5X,3.0X,3.5X |
| Battery type | Rechargeable lithium-ion (charger included) |
| Battery Life | 5-10HRS |
| Headlight housing | Aluminum |
| Headlight weight | 193g(With 1 pc Battery) |
| Headband material | ABS ratchet adjust ; antimicrobial protection on pads |
| Color Rendering Index | >93 |



LED SPOTLIGHT

Adjustable brightness for optimal visibility, ensuring versatility in a variety of applications.



LIGHTWEIGHT BATTERY

Lightweight and compact, it offers 5-12 hours of continuous power.



ADJUSTABLE HOOD

Elevate comfort for doctors with an optimized hood. Durable and adjustable, it caters to diverse needs.



COOL

Low intensity / Medium intensity / High intensity

ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales. The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

www.mediksurgical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US