


MK-H2400

LED Surgical Headlight



*Better care
Better comfort!*

MK-H2400

LED SURGICAL HEADLIGHT

DESCRIPTION

The MK-H2400 LED Surgical Headlight is specifically crafted for ENT, Dental Clinic, Vet, Urinary, Plastic surgery, Dermatology, and Hair Transplant, it offers unmatched brightness, adjustable spot size, and a durable 4400mAh battery. Ensure precision and clarity across diverse medical applications with this advanced lighting solution.

FEATURES

- The "PVC" Headband with aluminum material is safe and stable in quality.
- The latest Cold LED technology makes the illumination in better brightness max. 45,000Lux.
- Smoothly and uniform spot size in adjustable 4400MAH rechargeable battery
- Comparable with daylight in high color rendering index ≥ 93 .
- Clip-on which provides a longer working time of 5-10 hrs.
- The aluminum portable suitcase is better for protecting the headlight during shipment.

Model	MK-H2400
Luminance (LUX)	Up To 45,000
Illumination Spot diameter(D=Φ200MM)	10-78mm
Illumination Spot diameter(D=Φ300MM)	18-102mm
Illumination Spot diameter (D=Φ500MM)	35-180mm
Color temperature (K)	4000-5000
LED life	50000hrs
Battery Life	5-12 hrs
Color Rendering Index	>93
Magnification Glasses Optional	2.5X,3.0X,3.5X,4.0X,5.0X,6.0X
Battery type	Rechargeable lithium-ion (charger included)
Battery Life	5-12HRS
Headlight housing	Aluminum
Headlight weight	403g(With 1 pc Battery)
Headband material	ABS ratchet adjust; antimicrobial protection on pads
Color Rendering Index	>93



Adjustable Led Spotlight

Achieve optimal visibility and focus with adjustable spot size, ensuring versatility in various applications.



Lightweight Battery

Lightweight and compact, it offers 5-12 hours of continuous power.



Adjustable Pvc Hood

Elevate comfort for doctors with an optimized PVC hood. Durable and adjustable, it caters to diverse needs.



Cool

Low intensity / Medium intensity / High intensity .

ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales. The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

www.mediksurgical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US