

MK-H2600

Surgical Headlight



*Better care
Better comfort!*

MK-H2600

SURGICAL HEADLIGHT

DESCRIPTION

The MK-H2600 LED Surgical Headlight is tailor-made for a wide array of medical disciplines, including ENT, Dental Clinic, Veterinary, Urinary, Plastic Surgery, Dermatology, and Hair Transplantation. Featuring an integrated, head-mounted battery, it delivers unmatched brightness, adjustable spot size, and long-lasting performance, ensuring meticulous precision and clarity throughout various medical procedures.

FEATURES

- The "PVC" Headband with aluminum material is safe and stable in quality
- The latest Cold LED technology makes the illumination in better brightness max. 45,000Lux
- Comparable with daylight in high color rendering index ≥ 93
- Clip-on which provides a longer working time of 5-10 hrs
- Spot size diameter adjustable

Model	MK-H2600
Luminance (LUX)	Up To 45,000
Illumination Spot diameter(D=Φ200MM)	10-78mm
Illumination Spot diameter(D=Φ300MM)	18-102mm
Illumination Spot diameter (D=Φ500MM)	35-180mm
Color temperature (K)	4800
LED life	50000hrs
Battery Life	5-12 hrs
Color Rendering Index	>93
Magnification Glasses Optional	2.5X,3.0X,3.5X,4.0X,5.0X,6.0X
Battery type	Rechargeable lithium-ion (charger included)
Battery Life	5-12HRS
Headlight housing	Aluminum
Headlight weight	403g(With 1 pc Battery)
Headband material	ABS ratchet adjust; antimicrobial protection on pads
Color Rendering Index	>93



Adjustable Led Spotlight

Achieve optimal visibility and focus with adjustable spot size, ensuring versatility in various applications.



Multiple Adjustment Positions

The knob can be adjusted to different positions to fit the different angles required by the doctor.



Lightweight Battery

Embedded light battery, do not wear it on the body, the battery lasts long.



Adjustable PVC Hood

Elevate comfort for doctors with an optimized PVC hood. Durable and adjustable, it caters to diverse needs.



Cool

Low: intensity
 Medium: intensity
 High: intensity

ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales.The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

www.mediksurgical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US