

# MK-S45

Stainless steel instrument trolley



*Better care  
Better comfort!*

# MK-S45

# STAINLESS STEEL INSTRUMENT TROLLEY

## DESCRIPTION

MK-S45 multifunctional stainless steel instrument trolley is constructed with rectangular tube without any dead corner to accumulate the dust, easy for hygiene. The special fenced and reversed area is very practical for bottle or can container position. With functional accessories like storage drawer, dirt basin, trash can, etc., making this cart real practically for using, and can be used in first aid Room, general ward, ICU and other places.

## FEATURES

- 304 stainless steel material: stainless steel square tube or round tube, tube diameter 25-40mm, wall thickness 1.0-1.2mm, polished stainless steel countertop, thickness 0.8-1.0mm
- Casters: Silent wheels with a diameter of 100mm, two of which have counter brakes
- Contains stainless steel barrels, stainless steel trays, stainless steel basin brackets
- Contains 2 drawers for storing surgical instruments
- Guardrail on the countertop

## STANDARD SPECIFICATION

Size	90*50*95cm
Casters	100mm
Material	304 stainless steel

## ACCESSORIES

1 stainless steel clamp bracket	√
5 glass jars with lids	√
5 stainless steel tanks with lids	√
1 6-ounce stainless steel bowl	√
Stainless steel trash can	√



**Instrument tray**

Stainless steel movable instrument tray , It mainly stores scissors, tweezers and other instruments, which can be rotated and is convenient to use.



**Garbage can**

Rotatable garbage can with handle easy to carry ,stainless steel material is not easy to rust.



**Desktop disinfection sponge bottle holder**

There is a fixed shelf on the table top to prevent the bottle with the alcohol cotton from sliding.



**Push handle**

With handles, it is more convenient for nurses to push forward.



**Storage shelf & drawers**

The storage layer with guardrail can prevent the goods from sliding and falling during the movement, and there are 2 drawers to store things.



**Castor**

Silent wheels with counter brakes can be easily fixed after reaching the position.

## ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales. The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

**Pls scan the QR code for more info**

**[www.mediksurgical.com](http://www.mediksurgical.com) | +86 139 6223 1292**

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US