

MK-S46

Medical Instrument Trolley



*Better care
Better comfort!*

MK-S46

MEDICAL INSTRUMENT TROLLEY

DESCRIPTION

The MK-S46 Stainless Steel Hospital Instrument Trolley is a versatile and indispensable addition used in diverse health-care facilities. Engineered with precision and durability, this trolley is designed to streamline dressing procedures, enhance organization, and facilitate the mobility of medical supplies. With its thoughtful features including drawers, movable trays, and hooks, it offers a comprehensive solution for medical professionals to efficiently manage patient care tasks.

FEATURES

- Removable tray
- Stainless steel material is resistant to stains and easy to clean
- Four $\varnothing 75$ mm casters, two with brakes.
- Structure made of stainless steel
- Integrated Infusion bag hooks

STANDARD SPECIFICATION

Size	60*45*85cm
Material	Stainless steel
Castor	Dia75mm



Integrated Hooks

To further optimize utility, the trolley is outfitted with hooks for hanging infusion bags or instruments. This feature enhances accessibility and helps in maintaining a tidy workspace.



Movable Trays

The trolley features removable trays that are easily adjustable to accommodate various items. These trays can be moved and are easy to clean.



Drawer

The cart with the drawer to put the things ,it is designed for easy access and efficient organization.



Push Handle

With push handle on both sides which is easy to move the cart .

ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales.The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@medik-medical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US