

MK-E400ZL

Portable Examination Light



*Better care
Better comfort!*

MK-E400ZL

PORTABLE EXAMINATION LIGHT

DESCRIPTION

Elegant structure equipped with a multi-angle adjustable arm and a stable base, E400ZL is a mobile surgical lamp with OSRAM bulbs. Multiple adjustable parameters allow it to quickly adapt to various inspection and surgical scenarios.

FEATURES

- LED light source, almost no heat and long service life.
- Stable base and flexible arm, adjustable at multiple angles.
- Ultra-thin lamp head, easy to dissipate heat.
- Brightness, color temperature and spot diameter are adjustable.
- Removable sterilizable handle.
- Battery is optional.

TECHNICAL PARAMETER

Model	MK-E400ZL
Illuminance (Lux)	≥80000lux(adjustable)material,
Color temperature (K)	4000~5000K(adjustable)
Spot diameter (mm)	180~200mm(adjustable)
Color rendering index	Ra≥90
Bulb life	≥50000H



LED bulbs

Imported OSRAM lamp beads, with 50,000 hours of ultra-long lamp life, and equipped with ultra-thin lamp holders, excellent heat dissipation performance.



Sterilizable handle

Removable sterilizable handle for adjusting the spot diameter.



Modular design

Modular lighting system can reduce the fault area.



Function

Brightness, color temperature and spot diameter can be adjusted.



Flexible arm

Flexible joints can be adjusted at multiple angles.



Base

H-shaped base, strong and stable, double-sided silent castor with brakes.

ABOUT US

Medik is a group focusing on surgical room solution and nursing care solution

Medik is a group focusing on surgical room solution and nursing care solution. We started to manufacture surgical table & lamp since 2011 with a professional team consisting of experienced designer and engineer from international top brand. With 10 years accumulation, Medik is becoming a new star to provide reliable product and valuable solution for surgical room.

Qualified with ISO9001, ISO13485, CE, Medik has successfully distributed the surgical products covering surgical table & lamp & pendant to international hospitals both domestic and overseas including many countries in Europe.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@mediksurgical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US