

MK-E80JD

Examination Lamp



*Better care
Better comfort!*

MK-E80JD

EXAMINATION LAMP

DESCRIPTION

LED light source, compact design, easy to move and stay structure, all the excellent features make the MK-E80JL examination lamp suitable for most examination applications in hospitals, doctor's offices, and veterinary fields.

FEATURES

- Ideal for medical exams and procedures
- Adjustable arm levers allow for easy vertical and horizontal positioning
- Space-saving design folds to rest position against any wall
- Five, low friction, lockable casters on Floor Stand model ease positioning and securing placement
- Easy-to-clean smooth housing design aids infection prevention measures

STANDARD SPECIFICATION

Spot Diameter	15-270mm
LED Life	50000hrs
Color Temperature	4500K
Facula diameter	60-302mm
Light intensity	>35,000lux
Total Height	1800MM-1900MM
Working Distance	300MM-1400MM
LED Quantity	3pcs
Input	AC95-245V, 50/60Hz
Max Power	7W



LED lamp head

The examination lamp composed of 3 LED bulbs is suitable for most exam applications.



Gooseneck tube

Connected to the lamp head is a gooseneck, which can be bent at will.



Base

Five-grip aluminum alloy base, strong and stable(Optional plastic base).



Castor

Examination lamp with 5 castors, 2 of which have brakes.

ABOUT US

Medik is a group focusing on surgical room solution and nursing care solution. We started to manufacture surgical table & lamp since 2011 with a professional team consisting of experienced designer and engineer from international top brand. With 10 years accumulation, Medik is becoming a new star to provide reliable product and valuable solution for surgical room.

Qualified with ISO9001, ISO13485, CE, Medik has successfully distributed the surgical products covering surgical table & lamp & pendant to international hospitals both domestic and overseas including many countries in Europe.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@mediksurgical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US