

MK-E80JL1

Mobile Exam Light



*Better care
Better comfort!*

MK-E80JL1

MOBILE EXAM LIGHT

DESCRIPTION

The MK-E80JL1 surgical assistant illumination light is designed for optimal performance in the exam room. Featuring a uniform distribution of light spots, With 12W and 15W as optional, this versatile light source allows adjustable intensity. The flexible hose enables the adjustment of angles according to specific requirements, providing convenience and precision during examination.

FEATURES

- * Uniform light spot distribution and adjustable
- * Adjustable light intensity: 12W and 15W options
- * Flexible hose for customizable angle adjustments
- * Designed for enhanced surgical assistance

| | PARAMETER | PARAMETER |
|-----------------------|----------------|----------------|
| Power | 12W | 15W |
| Luminance (LUX) | 20000-50000lux | 20000-65000lux |
| Spot diameter(mm) | 15-270mm | 15-250mm |
| Color temperature (K) | 5000 ± 10% | 5000 ± 10% |
| LED life | 50000hrs | 50000hrs |



Adjustable Spot Distribution Control Knob

The intuitive knob interface enables easy adjustment of the spot distribution, providing a seamless user experience.



Foot-Operated Switch

The foot-operated master switch eliminates the need for surgeons to use their hands for light control, allowing them to maintain a sterile environment and focus on the surgical task.



Built-in Brightness Adjustment Panel

The brightness control offers variable intensity settings, allowing medical professionals to select the most suitable level of illumination for different stages of a procedure. This adaptability enhances the overall surgical experience.



Mobile Base

With a mobile base that allows easy movement and positioning of the illumination light.

ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales. The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

www.mediksurgical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US