

# MK-S01

Stainless steel Instrument trolley



*Better care  
Better comfort!*

# MK-S01

# STAINLESS STEEL INSTRUMENT TROLLEY

## DESCRIPTION

MK-S01 is a combined set instrument trolley which is designed for space saving due to the stack storage. Different height can be used in different environments and doctors.

## FEATURES

- 1.Nesting / stackable trolley
- 2.There are three trolleys in one set
- 3.The trolley is constructed with stainless steel material
- 4.Safety folds to all edges
- 5.One trolley with two shelves
- 6.With four 75mm silent castors and two with brakes

## STANDARD SPECIFICATION

Overall size (L*W*H)	1100*600*850/905mm900*550
	*800mm700*450*750mm
Castor	Dia 75mm



**Three different sizes**

There are three sizes in one cart , which will save a lot storage space . Also the different heights will satisfy the different needs from doctors .



**Side rail**

The side rail prevents the product from falling and can be used as a push handle.



**Polished technique**

Every corner has been carefully polished ,more rounded and beautiful.



**Castor**

The bedside table is equipped with dia 75 mm casters, 2 with brakes, for easy movement.

## ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales. The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

**Pls scan the QR code for more info**

**info@mediksurgical.com | +86 139 6223 1292**

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.