

# MK-S04

Stainless Steel Instrument Trolley



*Better care  
Better comfort!*

# MK-S04

# STAINLESS STEEL INSTRUMENT TROLLEY

## DESCRIPTION

MK-S04 is an instrument trolley with high-quality stainless steel, it ensures durability and resistance to corrosion, making it ideal for sterile environments. With its smooth-rolling wheels and sturdy construction, it guarantees reliable transportation of instruments, promoting seamless operations in busy hospital settings.

## FEATURES

- Smooth stainless steel shelves with handrail on three sides.
- Three stainless steel shelves with three sided railing on both shelf.
- Four 75mm silent castors and two with brakes are equipped.
- The guardrails are equipped to prevent objects from falling.
- Polished and welded well.
- Firm and stable installation.

## STANDARD SPECIFICATION

Material	Stainless steel
Shelves	2pcs
Castors	Four Dia 75mm castor two with brake
Size	900*550*800mm(Customize)



**Push Handle**

The stainless steel trolley features two handles on two sides, providing ease of maneuverability and effortless transportation.



**Polished Technique**

The polishing technique achieves a smooth angle and a flawless and lustrous finish, enhancing both the aesthetic appeal and durability of the material.



**Side Rail**

With 3 sides side rails to prevent things from falling while moving.



**Castor**

The stainless steel cart is equipped with dia 75 mm casters, 2 with brakes, for easy movement.