

MK-S17B

Stainless Steel linen Cart



*Better care
Better comfort!*

MK-S17B

STAINLESS STEEL LINEN CART

DESCRIPTION

MK-S17B is a linen cart with the high quality stainless steel ,it is mainly used in hospitals for finishing bed sheets and pillowcases.It is greatly improved the work efficiency of nurses.

FEATURES

- Push handle at one end.
- Stainless steel frame
- Big bag for wasterNoiseless operation
- Four Dia75mm casters and two with brake

STANDARD SPECIFICATION

Overall size (L*W*H)	500*450*950mm
Castor	Dia 75mm



Plastic Top Frame

Plastic top frame is easier to clean and more aesthetically pleasing.



Push Handle

The cart with push handle, more convenient for medical staff to operate.



Bag Material

Made of Oxford cloth, detachable and easy to wash.



Castor

The stainless steel cart is equipped with dia 75 mm casters, 2 with brakes, for easy movement.