

# MK-S30B

Stainless Steel  
Medical Kick Bucket



*Better care  
Better comfort!*

# MK-S30B

# STAINLESS STEEL MEDICAL KICK BUCKET

---

## DESCRIPTION

Heavy-duty construction, thick materials, non-slip details, and stable casters allow the MK-S30B kick bucket to operate smoothly during surgery, and the anti-tip, rubber chassis protects the wall.

---


## FEATURES

- High grade heavy duty stainless steel structure
- Anti-slip rubber bumper ring and anti-slip fixing points
- Four 50mm smooth castors with brakes
- Anti-dumping design, smooth operation

---

## STANDARD SPECIFICATION

Safe Working Load	15L
Castor	Dia 50mm
Material	Stainless steel



**Rubber bumper ring**

Anti-slip rubber bumper ring to protect the walls.



**Anti-slip fixing points**

2 units of rubber padded on each column to fix the bucket.



**Bucket**

Removable full stainless steel bucket ,heavy duty .



**Castors**

4 units silent castors with individual brakes.

## ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales.The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

**Pls scan the QR code for more info**

**info@medik-medical.com | +86 139 6223 1292**

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.