

MK-33B/MK-34B

Hand Wash Basin Trolley



*Better care
Better comfort!*

MK-33B MK-34B

HAND WASH BASIN TROLLEY

DESCRIPTION

MK-33B and MK-34B is a Basin Stainless Steel Hand Basin trolley, this innovative solution combines functionality, hygiene, and convenience, all in one sleek package. It can be used in different hospital departments with dual basins, an integrated shelf and castors.

FEATURES

- Hand washing/storage trolley
- The trolley is constructed with stainless steel material
- Safety folds to all edges
- One trolley with one shelf and wash basin
- With four 50 mm silent castors and two with brakes

STANDARD SPECIFICATION

Model	MK-S33B	MK-S34B
Overall size (L*W*H)	380*800*850mm	380*400*850mm
Castor	Dia 50mm	Dia 50mm
Basin capacity	10L	5 L

**Wash Basin**

The 5L stainless steel wash basin is durable, hygienic, easy to clean, ample capacity and corrosion resistant.

**Storage Shelf**

The trolley with the shelf to storage the items ,It increases storage capacity and saves space.

**Polished Technique**

Every corner has been carefully polished ,more rounded and beautiful.

**Castor**

The stainless steel cart is equipped with dia 75 mm casters, 2 with brakes, for easy movement.

ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales.The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@medik-medical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US