

MK-S49

Instrument Trolley



*Better care
Better comfort!*

MK-S49

INSTRUMENT TROLLEY

DESCRIPTION

Square tube reinforcement structure at the bottom of the trolley, thick shelf and perfect polished finishing, flexible castors, make MK-S49 can be applied to a variety of scenarios. Also the knock down design can save cost.

FEATURES

- Two stainless steel shelves.
- Smooth stainless steel shelves with handrail on two sides.
- Four 100mm silent castors and two with brakes are equipped.
- Polished and welded well.
- Additional square tube reinforcement structure basement, stable and durable.

STANDARD SPECIFICATION

Overall Size	L810*W450*H900mm
Safe Working Load	100kg
Castor	Dia 100mm
Packing	Knock down design
Structure	Full stainless steel




Push handle

The stainless steel trolley features two handles on two sides, providing ease of maneuverability and effortless transportation.



Polished technique

The polishing technique achieves a smooth angle and a flawless and lustrous finish, enhancing both the aesthetic appeal and durability of the material.



Basement

Square tub reinforcement structure at the bottom of the trolley will make it durable and stable.



Castors

4 units silent castors with diagonal lock, easy to move.

ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales. The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

www.mediksurgical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US