

MK-S50

Instrument Trolley



*Better care
Better comfort!*

MK-S50

INSTRUMENT TROLLEY

DESCRIPTION


Perfectly polished and welding, smooth interfaces and edges, rich configurations and sufficient storage space, full stainless steel material,antibacterial,sterilizable and cleanable, make this instrument cart suit for various usage scenarios in the hospital.

FEATURES


- Full stainless steel structure ,antibacterial and easy to clean
- Three stainless steel shelves with guardrail.
- Four castors with bumpers and diagonal lock
- Perfectly polished finishing
- With adjustable tray ,bucket,bowl and drawer
- With bottle holders

STANDARD SPECIFICATION

Overall Size	L 770*W 470*H 800mm
Tray size	L 260*W 180mm
Bucket size	10 L
Packing	Knock down design
Structure	Full stainless steel
Castor	Dia 75mm



Bottle holder on the table top.



Adjustable stainless steel bowl.




Adjustable stainless steel tray.



Adjustable stainless steel bucket.



Stainless steel drawer for storage.



4 units silent castors with bumpers and diagonal lock , easy to move.

ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales.The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@medik-medical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US