

YA-01E

Surgical table



*Better care
Better comfort!*

YA-01E

SURGICAL TABLE

DESCRIPTION

Easy transportation and re-location, superior stability due to unique structure, advanced position parameters, excellent finishing details, all these features make the surgical table YA-01E as an ideal option for doctors for various surgery like thoracic surgery, abdominal surgery, gynecology and obstetrics, ophthalmology and ENT, urology, colon surgery.

MAIN FEATURES

- Mechanical system for waist board lifting of electric operating table
- Gas spring for adjusting movements of headrest section
- Polyurethane shaped mattress which is easy to remove and clean
- Four 125mm wheels for moving the OT table, which are with central braking system
- AISI304 stainless steel structure and chassis, which is easy to clean
- Polyurethane shaped mattress which is easy to remove and clean
- With movements of horizontal and longitudinal for C-arm photography examination
- With four 125mm wheels to move the table, which are with twin-wheel central braking system
- Eight independent motors for movements of HI-LO, back section, Trendelenburg and Reverse-Trendelenburg, Lateral tilting function, Horizontal and longitudinal movement (For C-arm examination), Leg section (separated left or right leg section, also left and right leg section as one united set), auto-contour of operating table
- Gas spring for headrest movement
- Mechanism system for waist board lifting
- Leg section can be removable, then connected to orthopedic traction shelf directly
- With lateral control system incorporated to side column, which is with blocking system for movements

Main details



Remote controller

Remote control can backlit display indicating status information, it ensures a safe and easy handling even in the dark operation room.



Column controller

The column controller is also with backlit, it is more convenient and fast to operate.



Kidney bridge

The kidney bridge is adjusted by the crank which can up to 120mm.



Form mattress

High-quality poly surface for easy clean, removable, impermeable



Stainless steel basement

The basement is made of stainless steel 304, which is stable and with a good outlook.



5" big Castor

5" wheels make it easy to move the surgical table to different locations.

Technical parameters

Model	YA-01E
Table Length	2200mm (88")
Table Width	600mm(24")
Table Height Range	700-1000mm (28" to 40")
Kidney elevator mechanism	120mm(4.8")
Slide motion	350mm(14")
Back section up/down	70°/22°
Head section up/down	45°/90°
Leg section up/down	15°/90°
Reverse trendelenburg	30°
Trendelenburg	30°
Lateral tilt left/right	22°/22°
Load capacity	320kg

Standard accessory

Headrest	1pc
Armrest	1 pair
Waist support holder	1 pair
Shoulder support holder	1 pair
Leg support holder	1 pair
Footrest	1 pair
Anesthesia screen frame	1 pc
Body strap	1 set
Handset controller	1 pc
Lateral controller inner constructed on side of table	1 set

Professional functions

Lowest height position

700MM



Highest height position

1000MM



Back section motion

-22° to 70°



Slide motion

350MM



Reverse trendelenburg

30°



Trendelenburg

30°



Lateral tilt left

22°



Lateral tilt right

22°



Professional functions

Head section up

45°



Head section down

90°



Leg section motion

-90° to +15°



Leg section outward

90°



Kidney elevator mechanism

120MM



ABOUT US

Medik is a group focusing on surgical room solution and nursing care solution. We started to manufacture surgical table & lamp since 2011 with a professional team consisting of experienced designer and engineer from international top brand. With 10 years accumulation, Medik is becoming a new star to provide reliable product and valuable solution for surgical room.

Qualified with ISO9001, ISO13485, CE, Medik has successfully distributed the surgical products covering surgical table & lamp & pendant to international hospitals both domestic and overseas including many countries in Europe.



YOU TUBE



MEDIK

Pls scan the QR code for more info

www.mediksurgical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US