

YA-05E

Surgical Table



*Better care
Better comfort!*

YA-05E

SURGICAL TABLE

Engineered for increasingly complex surgical scenes, YA-05E comes with intelligent automatic control system satisfying various demands of general surgeries. Extra-wide tabletop, extreme loading system, long horizontal sliding that can be suitable for both X-ray and C-arm. YA-05E is an ideal electric operating table.

FEATURES

- Four 125mm wheels with centralized brake system
- Divided 6 sections top plate with various patient positions for all kinds of general surgery
- Horizontal table top super 350mm sliding for C-arm access to entire body
- With access to connect the extension device and traction for orthopedic operations
- Double deck table top for X-ray cassette position
- With a button to operate the table in one step, return the table to a completely horizontal or "zero" position
- With memory function of the movements for the doctor to find the position for certain surgery operation in short time
- Built in kidney bridge
- Memory mattress with ultrasonic seamless welding cover to provide great comfort for long time surgery
- Super SS frame and base anti corrosion
- Hand remote control with LED light integrated to the manual control to illuminate it when the operating room is darkened
- Inbuilt column control panel with visible icon for position operation

ITEM	PARAMETER
Table Length	2060mm
Table Width	550mm
Height adjustment	700-1000mm
Head plate adjustment	+40°to -90
Back plate adjustment	+80°to -20°
Leg plate adjustment	+15°to -90°
Leg plate outward	90°
Trendelenburg	25°
Reverse trendelenburg	25°
Lateral Tilt (left/right)	15°/ 15°
Kidney bridge	110mm
Sliding	350mm
Loading capacity	250kg

Lowest height position
700mm



Highest height position
1000mm



Back plate up
0-80°



Back plate down
0-20°



Reverse trendelenburg
0-25°



Trendelenburg
0-25°



Lateral tilt left
0-15°



Lateral tilt right
0-15°



Horizontal sliding
350mm



Head section up
0-40°



Head section down
0-90°



Leg section down
0-90°




Leg section outward
0-90°



Kidney elevator
110mm





Material

The base material is stainless steel304, which is easier to clean, modern and stable



Power

internal battery maintaining power supply for 50-80 operations



kidney bridg

built-in kidney bridge, more convenient




5" castors central brake




Handle controller and column controller



Body support with cushion



Leg support with cushion



Arm Support with Cushion



Anesthesia Screen Frame



Shoulder Support with Cushion

ABOUT US

Medik is a group focusing on surgical room solution and nursing care solution

Medik is a group focusing on surgical room solution and nursing care solution. We started to manufacture surgical table & lamp since 2011 with a professional team consisting of experienced designer and engineer from international top brand. With 10 years accumulation, Medik is becoming a new star to provide reliable product and valuable solution for surgical room.

Qualified with ISO9001, ISO13485, CE, Medik has successfully distributed the surgical products covering surgical table & lamp & pendant to international hospitals both domestic and overseas including many countries in Europe.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@mediksurgical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US