

YA-10EB

Surgical table



*Better care
Better comfort!*

YA-10EB

SURGICAL TABLE

YA-10EB is an electric surgical table which is suitable for the major universal surgeries like thoracic surgery, abdominal surgery, brain surgery, ophthalmology, otolaryngology, obstetrics and gynecology, urology, orthopedics, etc. Stainless steel material and heavy duty structure makes the table durable in use for long lifespan. The table with horizontal sliding function is suitable for specific surgery which needs whole X-ray translucent compatible for C-arm exam.

FEATURES

- Eight independent motors for movements of hi-lo, back section, trendelenburg and reverse-trendelenburg, lateral tilting function, horizontal and longitudinal movement(for c-arm examination), leg section(separated left or right leg section, also left and right leg section as one united set), auto-contour of operating table
- The table can be moved vertically, and it can be used with c-arm for radiographic diagnosis or filming.
- The countertop is made of high-strength melamine resin sheet, and the mattress is made of memory foam mattress.
- The leg plate can be removed, manually rotated abduction, folded down, easy to adjust, very convenient for urological surgery.
- The handheld remote control which is easy to operate, safe and reliable.
- Gas spring for adjusting movements of headrest section
- Polyurethane shaped mattress which is easy to remove and clean
- Aisi304 stainless steel structure and chassis, which is easy to clean
- Gas spring for headrest movement
- Mechanism system for waist board lifting
- Leg section can be removable, then connected to orthopedic traction shelf directly

ITEM	PARAMETER
Table Length	2010mm(79")
Table Width	500mm(20")
Table Height Range	700~1000mm(27.5"~39")
Back section up/down	75°/10°
Head section up/down	50°/90°
Leg section down/outward	90°
Reverse trendelenburg	25°
Trendelenburg	25°
Lateral tilt left/right	20°
Slide motion	350mm (13.5")
Kidney elevator mechanism	110mm (4.3")
Load capacity	250kg (550lb)

ITEM	PARAMETER
Headrest	1pc
Armrest	1 pair
Waist support holder	1 pair
Shoulder support holder	1 pair
Leg support holder	1 pair
Footrest	1 pair
Anesthesia screen frame	1 pc
Body strap	1 set
Handset controller	1 pc
Lateral controller inner	1 set

constructed on side of table

Lowest height position
700mm



Highest height position
1000mm



Back plate up
0-75°



Back plate down
0-10°



Reverse trendelenburg
0-25°



Trendelenburg
0-25°



Lateral tilt left
0-20°



Lateral tilt right
0-20°



Horizontal sliding
350mm



Head section up
0-50°



Head section down
0-90°



Leg section down
0-90°



Leg section outward
0-90°



Kidney elevator
110mm





Remote controller

The remote control icon is clear and large, which is convenient for quick identification and operation.



Kidney bridge

The kidney bridge is adjusted by the crank which can up to 110mm



Form mattress

High-quality poly surface for easy clean, removable, impermeable



Stainless steel basement

The basement is made of stainless steel 304, which is stable and with a good outlook.



Castors and lock system

The surgical table is equipped with casters easy to move the surgical table to the place they need , and it with the lock system.

ABOUT US

Medik is a group focusing on surgical room solution and nursing care solution

Medik is a group focusing on surgical room solution and nursing care solution. We started to manufacture surgical table & lamp since 2011 with a professional team consisting of experienced designer and engineer from international top brand. With 10 years accumulation, Medik is becoming a new star to provide reliable product and valuable solution for surgical room.

Qualified with ISO9001, ISO13485, CE, Medik has successfully distributed the surgical products covering surgical table & lamp & pendant to international hospitals both domestic and overseas including many countries in Europe.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@mediksurgical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US