

YA-11E

Surgical table



*Better care
Better comfort!*

YA-11E / SURGICAL TABLE

Extreme flatten base, easy to clean whole frame, YA-11E with the unique features is a new generation surgical table designed to have a better stability, more access for surgeon to patient, wider space for C-arm operation, etc..

Main Features

- The contact area between the base and the ground is large and the stability is extremely strong.
- AISI304 stainless steel construction and chassis, easy to clean
- Polyurethane-shaped mattress, easy to remove and clean
- With horizontal and vertical movement, used for C-arm photography inspection
- Lateral control system integrated into the side column with motion blocking system
- Gas spring for adjusting movements of headrest section
- Equipped with a trolley to move the operating table
- Eight independent motors for movements of HI-LO, back section, Trendelenburg and Reverse-Trendelenburg, Lateral tilting function, Horizontal and longitudinal movement(For C-arm examination), Leg section(separated left or right leg section, also left and right leg section as one united set), auto-contour of operating table
- Mechanism system for waist board lifting

Technical parameters

| | |
|---------------------------|-----------------------|
| Model | YA-11E |
| Table Length | 2150mm(84.5") |
| Table Width | 550mm(21.5") |
| Table Height Range | 720-1070mm(28"~42") |
| Back section up/down | 72°/32° |
| Head section up/down | 30°/90° |
| Leg section down | 90° |
| Reverse trendelenburg | 25° |
| Trendelenburg | 25° |
| Lateral tilt left/right | 22° |
| Slide motion | 300mm (12") |
| Kidney elevator mechanism | 120mm (5") |
| Load capacity | 275kg (600lb) |

Standard accessory

| | |
|----------------------|----------|
| Mattress | 1 set |
| Anesthesia screen | 1 piece |
| Shoulder support | 1 pair |
| Arm board | 1 pair |
| Body support | 1 pair |
| Leg support | 1 pair |
| Foot plate | 1 pair |
| Kidney bridge handle | 1 piece |
| Hand remote | 1 piece |
| Fixing clamp | 8 pieces |

Main details



Basement

The base is in contact with the ground over a large area and has an extreme stability. Meanwhile, the surgeon get more access to the patient with flat base. Moreover, the movement space of the C-arm is increased, which increases the operation space of the C-arm specially for those big machine.



Remote controller

Remote control can backlit display indicating status information, it ensures a safe and easy handling even in the dark operation room.



Column controller

The column controller is also with backlit, It is more convenient and fast to operat.



Mattress

Polyurethane molded mattress, waterproof, easy to disinfect and clean. Snap-fastened and easy to remove.



Mobile Trolley

The mobile trolley is customized to lift and mobilize the table, with minimum operation, it can lift the table onto the trolley and transfer to a wanted place.

Professional functions

Lowest height position
720mm



Highest height position
1070mm



Back section up motion
72°



Back section down motion
32°



Head section up
30°



Head section down
90°



Leg section motion
-90°



Leg section outward
90°



Slide motion
350mm



Kidney elevator mechanism
120mm



Reverse trendelenburg
30°



Trendelenburg
30°



Lateral tilt left
22°



Lateral tilt right
22°



ABOUT US

Medik is a group focusing on surgical room solution and nursing care solution. We started to manufacture surgical table & lamp since 2011 with a professional team consisting of experienced designer and engineer from international top brand. With 10 years accumulation, Medik is becoming a new star to provide reliable product and valuable solution for surgical room.

Qualified with ISO9001, ISO13485, CE, Medik has successfully distributed the surgical products covering surgical table & lamp & pendant to international hospitals both domestic and overseas including many countries in Europe.



YOU TUBE



MEDIK

Pls scan the QR code for more info

INFO@MEDIKSURGICAL.COM | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US