

# YA-B103

## Orthopedic Traction Frame



*Better care  
Better comfort!*

# YA-B103

# ORTHOPEDIC TRACTION FRAME

## DESCRIPTION

The YA-B103 orthopedic traction extension with carbon fiber and double joint design can play the best performance in perspective and adjustment range. It is suitable for C-arm and can provide doctors with the maximum visual range. The multi-dimensional and multi-angle adjustment of each part can assist doctors to obtain the best surgical position. The carbon fiber material can let the C-arm generate a perfect fluoroscopy image during surgery.

## FEATURES

- Main part is made of carbon fiber, which has strong perspective and is suitable for X-ray and C-arm.
- Dual joint design ,larger adjustment range and wider surgical field of view.
- The boot can be adjusted in various heights, angles and distances to provide the best surgical position.
- The traction mechanism has a precise scale display.
- Equipped with an auxiliary trolley, it is easy to install and use, and is suitable for various operating tables.

## STANDARD SPECIFICATION

Height adjustment	Can follow the height of surgical table
Carbon Fiber Pole Outreach	360°
Traction stroke	≥180mm
Horizontal movement range of knee support	≥550mm
Vertical angle adjustment range of traction mechanism	50°~140°
Overall size	1800*350mm (Width and height can be adjusted to the operating table)



**Carbon fiber material**

The seat board and the main part are made of carbon fiber material, which is suit for X-ray and C-arm, no blind spots in any surgical position



**Auxiliary trolley**

Equipped with an auxiliary trolley for easy installation and use



**Knee support**

Equipped with a knee support to increase flexibility during surgery



**Foot Adjustment Device**

The foot section can be adjusted in multiple angles and distance ,as well as the height ,allowing doctors to easily fix the patient in the best surgical position



**Traction crank**

The traction crank can be used to traction and stretch the joints after fixing the patient's foot, With precise scale indication,making it easier for doctors to perform surgery



**Height adjustment device**

The height of the traction frame can be adjusted according to the height of the operating table so that they remain balanced



**Dual Joint**



Double-joint traction frame, with a wider adjustable range, enables more surgical positions and provides the doctor with the largest surgical field of view

## ABOUT US

Medik is a group focusing on surgical room solution and nursing care solution

Medik is a group focusing on surgical room solution and nursing care solution. We started to manufacture surgical table & lamp since 2011 with a professional team consisting of experienced designer and engineer from international top brand. With 10 years accumulation, Medik is becoming a new star to provide reliable product and valuable solution for surgical room.

Qualified with ISO9001, ISO13485, CE, Medik has successfully distributed the surgical products covering surgical table & lamp & pendant to international hospitals both domestic and overseas including many countries in Europe.



YOU TUBE



MEDIK

Pls scan the QR code for more info

**info@mediksurgical.com | +86 139 6223 1292**

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US