

# YA-S09

Electric Doctor Stool For Surgery



*Better care  
Better comfort!*

# YA-S09

## ELECTRIC DOCTOR STOOL FOR SURGERY

### DESCRIPTION

The powered column lifting, ergonomically designed backrest and armrests, flexible adjusting arm support, the innovative surgical stool can provide doctors with stable and comfortable support during ophthalmology operations and medical examinations, making it easier for doctors to find their focus, improve work efficiency.

### MAIN FEATURES

- Ergonomic design for easy use
- Stable electric column structure
- Removable backrest with adjustable angle
- The armrests are height-adjustable and can be rotated 180 degrees
- Height of the chair can be adjusted easily and conveniently through the foot pedal
- Built-in battery, rechargeable
- Four castors for easy transport

### TECHNICAL PARAMETER

Model	YA-S09
Height range	460-600mm
Seat top diameter	Dia 420mm
Armrest adjustment	rotated 180 degrees
Castor	Dia 50mm

## PRODUCT DETAILS

### Back rest adjustment

Back rest can be adjusted by lifting handle and can be removed through the knob.



### Arm rest adjustment

Height of the arm rest can be adjusted by button, also can be rotated 180 degrees.



### Height adjustment

Height of the chair is adjusted by foot pedal.



### Charging interface


There is a charging interface and switch on the back of the chair base, and a built-in battery.



## PRODUCT DETAILS

### Column structure

Column structure, ensure stability during surgery.



---

### Base and Castors

4 swivel castors and stable base ,durable and easy to move.



## ABOUT US

Zhangjiagang Medi Medical Equipment Co.,Ltd, founded in 1995, with more than 20 years experience in developing hospital furniture, is one of the largest manufacturer of hospital bed, patient stretcher, medical cart and other related hospital furniture in China or even Asia.

Through scientific management complying with ISO9001&ISO13485, by adopting advanced facility like Germany TRUMPF CNC, Japan OTC Robot Welding, Auto-power coating line with Swiss GEMA gun , Medi focuses on providing brand quality product with CE approval.

Headquartered in Zhangjiagang nearby Shanghai, with subsidiary oversea offices in Hungary, Dubai, Thailand, New Zealand, Kenya, Colombia, Medi acting professionally and swiftly for both pre-sales and after-sales. The products are widely accepted by global clients, such as UK RIE hospital, UNIVERZITETNI KLINICNI CENTER, Kokilaben Dhirubhai Ambani Hospital, Singapore National University Hospital etc.



YOU TUBE



MEDIK

**Pls scan the QR code for more info**

**[www.mediksurgical.com](http://www.mediksurgical.com) | +86 139 6223 1292**

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US