

YA-SR2D

Double Station Wash Sink



*Better care
Better comfort!*

YA-SR2D

DOUBLE STATION WASH SINK

DESCRIPTION

YA-SR2D is a floor-standing double station sink, sink body is made of acrylic. It's unique curved appearance and corrosion-resistant materials make it suitable for medical staff's usage habits while being durable and splash-proof.

FEATURES

- Elegant appearance and easy to clean.
- The sink body is made of high-density acrylic material, which is corrosion-resistant, heat-resistant, and not easy to crack.
- Floor-standing curved sink, suitable for medical staff's usage habits.
- Inner stainless steel structure, heavy duty.
- Optional induction soap box, integrated LED light mirror, mechanical or electronic constant temperature water mixing regulating valve.

STANDARD SPECIFICATION

Overall Size L1500*H1230*D680mm
Sink High-density non-porous material, one-piece hot pressing
Structure Stainless steel



Integrated LED light mirror

Convenient lighting, high-definition mirror.



Sensor faucet

Induction water, water temperature can be adjusted freely.

ABOUT US

Medik is a group focusing on surgical room solution and nursing care solution

Medik is a group focusing on surgical room solution and nursing care solution. We started to manufacture surgical table & lamp since 2011 with a professional team consisting of experienced designer and engineer from international top brand. With 10 years accumulation, Medik is becoming a new star to provide reliable product and valuable solution for surgical room.

Qualified with ISO9001, ISO13485, CE, Medik has successfully distributed the surgical products covering surgical table & lamp & pendant to international hospitals both domestic and overseas including many countries in Europe.



YOU TUBE



MEDIK

Pls scan the QR code for more info

info@mediksurgical.com | +86 139 6223 1292

Medik reserves the right to make changes without notice in design, specifications and models. The only warranty Medik makes is the express written warranty extended on the sale or rental of its products.

© 2022 Medik Services, Inc. ALL RIGHTS RESERVED. APR53703 rev 1 01-JAN-2022 ENG – US